



# Glencar

## Activity Chair by Rifton



### TECHNICAL SPECIFICATIONS

	Small Standard Base	Small Hi/Lo Base	Medium Standard Base	Medium Hi/Lo Base	Large Standard Base	Large Hi/Lo Base
Frame Width	Short legs: 21	26	Short legs: 23 <sup>1</sup> / <sub>4</sub>	27 <sup>1</sup> / <sub>2</sub>	Short legs: 26 <sup>3</sup> / <sub>4</sub>	29 <sup>3</sup> / <sub>4</sub>
	Long legs: 23		Long legs: 24 <sup>1</sup> / <sub>2</sub>		Long legs: 26 <sup>3</sup> / <sub>4</sub>	
	Short legs with casters: 22 <sup>1</sup> / <sub>2</sub>		Short legs with casters: 24 <sup>1</sup> / <sub>2</sub>		Short legs with casters: 27 <sup>3</sup> / <sub>4</sub>	
	Long legs with casters: 23 <sup>1</sup> / <sub>2</sub>		Long legs with casters: 25 <sup>1</sup> / <sub>2</sub>		Long legs with casters: 27 <sup>3</sup> / <sub>4</sub>	
Seat Height Above Floor	Short legs: 9 <sup>1</sup> / <sub>2</sub> – 12 <sup>1</sup> / <sub>2</sub>	10 – 25	Short legs: 12 <sup>1</sup> / <sub>2</sub> – 16 <sup>1</sup> / <sub>2</sub>	12 – 23	Short legs: 16 <sup>1</sup> / <sub>2</sub> – 20 <sup>1</sup> / <sub>2</sub>	16 – 24
	Long legs: 18 <sup>1</sup> / <sub>2</sub> – 21 <sup>1</sup> / <sub>2</sub>		Long legs: 19 – 23		Long legs: 20 – 24	
	Short legs with casters: 13 <sup>1</sup> / <sub>2</sub> – 16 <sup>1</sup> / <sub>2</sub>		Short legs with casters: 16 – 20		Short legs with casters: 16 <sup>1</sup> / <sub>2</sub> – 20 <sup>1</sup> / <sub>2</sub>	
	Long legs with casters: 18 <sup>1</sup> / <sub>2</sub> – 21 <sup>1</sup> / <sub>2</sub>		Long legs with casters: 19 – 23		Long legs with casters: 20 – 24	
Seat Angle for Tilt-in-space	15° forward, 15° back	15° forward, 25° back	15° forward, 15° back	15° forward, 25° back	15° forward, 15° back	15° forward, 25° back



Glencar [Enquire about this product](#)

T: 01 8900201 or 028 41754854 E: [info@glencar.ie](mailto:info@glencar.ie) [www.glencar.ie](http://www.glencar.ie)

*grow, play, thrive*



# Glencar

## Activity Chair by Rifton



### TECHNICAL SPECIFICATIONS

	Small Standard Base	Small Hi/Lo Base	Medium Standard Base	Medium Hi/Lo Base	Large Standard Base	Large Hi/Lo Base
Backrest Angle Adjustment	5° forward, 20° back	5° forward, 20° back	5° forward, 20° back	5° forward, 20° back	5° forward, 20° back	5° forward, 20° back
Footboard Knee Angle	45° – 110°	45° – 110°	45° – 110°	45° – 110°	45° – 110°	45° – 110°
Seat Height Above Footboard (in)	7 – 12	7 – 12	7¾ – 15¼	7¾ – 15¼	11¼ – 21¼	11¼ – 21¼
Seat Width with Hip Guides (in)	7 – 9	7 – 9	8½ – 11½	8½ – 11½	11 – 14	11 – 14
Seat Width without Hip Guides (in)	12	12	14	14	17	17
Seat Depth (1" increments)	8 – 12	8 – 12	11 – 16	11 – 16	15 – 20	15 – 20
	7 – 11 with mini kit	7 – 11 with mini kit				
Armrest Height Above Seat (in)	5 – 7½	5 – 7½	7 – 10½	7 – 10½	7 – 10½	7 – 10½
Trunk Support Width (in)	5½ – 11½	5½ – 11½	6½ – 12	6½ – 12	9½ – 14	9½ – 14
Backrest Height Above Seat (in)	12½ – 15½	12½ – 15½	15½ – 19	15½ – 19	19 – 24	19 – 24
Headrest Height Above Heat (in)	14½ – 21	14½ – 21	17½ – 24½	17½ – 24½	19½ – 29½	19½ – 29½
Maximum Working Load (lbs)	75	75	150	150	250	250



Glencar [Enquire about this product](#)

T: 01 8900201 or 028 41754854 E: [info@glencar.ie](mailto:info@glencar.ie) [www.glencar.ie](http://www.glencar.ie)

*grow, play, thrive*